



# Leica HCS A

Amplify the Power of Imaging  
High Content Screening Automation


Living up to Life

*Leica*  
MICROSYSTEMS



- Amplify the power of imaging with Leica HCS A
- Easy-to-use automation provides efficient high content screening
- Maximum flexibility for universal applications
- Powerful hardware for high performance imaging





High Content Screening (HCS) allows researchers to quickly change from qualitative to quantitative fluorescence imaging during an experiment. Automated high resolution imaging therefore answers complex questions in less time. It simplifies research work and efficiently reveals relationships within and between cells and organisms. Leica Microsystems offers a set of innovative tools to convert your high resolution microscope into a high content imaging device.

## Leica HCS A


### High Content Screening Automation

Leica Microsystems provides a wide range of confocal and widefield systems, known for brilliant image quality and maximum performance. Combining high content screening with intelligent microscope automation greatly amplifies the power of an imaging system.

The value of an imaging system becomes more than the sum of its parts when the MATRIX M3 screening software is added to the Leica LAS AF platform (Leica Application Suite Advanced Fluorescence). More experiments can be performed by automated sample screening resulting in standardized experiment results. Quantification easily provides statistically relevant results.

Leica Microsystems' automated high content screening speeds up experiment throughput and enhances laboratory capacity. Automation reduces routine microscopy and improves workflow. From automated routine image acquisition to complex HCS experiments with on-the-fly image analysis, Leica HCS A is the right solution. Leica HCS A fully integrates with your laboratory environment via interfaces to existing image analysis packages. The Computer Aided Microscopy (CAM) tool enables free programming of the imaging system and the creation of specific protocols and workflows.

Customize the Leica microscope according to your actual needs and discover the unrivaled application flexibility for system biology, cancer research or environmental screening.




Modern research is a continuous cycle of experiment design, data acquisition and data handling, to answer questions about life's processes.



### Features

- High performance imaging
- Time saving automation
- Open architecture
- Platform independent results
- OME data formats
- Perfect integration

## Amplify the Power of Imaging




### Automated Leica High Content Screening



High resolution imaging techniques answer many questions in modern life science. For high content screening, automation is essential if researchers are to efficiently achieve results.

### Intelligent Automation

Leica HCS A adds extensive automation capability to confocal and widefield microscopes and converts stand alone systems into fully-featured high content screening devices.

MultiPosition-MultiParameter experiment designs with autofocus and drift compensation provide maximum imaging flexibility to obtain concise imaging results. The system benefits from open interfaces for full integration into your imaging facility. Leica generates Open Microscopy Environment (OME) and standard TIFF data for platform independent analysis. Existing algorithms and image analysis programs from open source or shared resources can be utilized, saving costs and time. External programming languages from all operating systems can address the acquisition system, giving flexible control for more demanding tasks.






**a** Rat brain slice, small neuron network layer 5. Interneurons (Alexa 594, red) and Pyramidal Cell Oregon (Bapta 1, calcium sensitive, green). Courtesy of Dr. Thomas Nevian, Institute of Physiology, University of Bern, Switzerland.

**b** *Danio rerio* – Zebrafish – Nuclear and Acetylated  $\alpha$ -Tubulin staining of 6 days flh:eGFP Zebrafish larvae Nuclei (Hoechst, blue), acetylated tubulin (red) and neurons (GFP, green). Courtesy of ICI Imaging Centre IGBMC, Illkirch, France.

For full flexibility in your experimental approach you can choose from selecting a predefined template or creating your own sophisticated protocol.

Since the fully automated experiment can be applied to a massive number of samples, acquisition speed is dramatically improved. The open architecture of Leica HCS A provides a direct link between the ongoing acquisition process and the image analysis, allowing for changes in imaging parameters such as the detection of a rare event. Acquisition speed, intelligent automation plus seamless laboratory integration maximize the power of Leica Microsystems' imaging systems and your research.




Many companies offer dedicated imaging routines, but only for dedicated assays. Leica Microsystems provides standard solutions for routine experiments plus maximum flexibility of design to suit your requirements.

# Easy-to-use Automation

## We Keep It Simple!

Wizards guide the user through an experiment in a streamlined way. Design follows function – benefit from clear user interfaces, ensuring fast training and the highest productivity.



### Predefined Scanning Templates

Place the specimen carrier on the microscope stage, enter the experiment ID, and move to the start point. Upload a pre-configured scanning template and fine-tune the screening job according to the experiment needs. With a single click of a button, all positions are automatically calculated and the experiment is ready to start. The image data is continuously streamed to a local or network attached storage device, ready for immediate analysis.


### Perfectly Timed Workflow

The operator stays in full control of the experiment and may stop or pause the screening at any time. The system provides continuous user feedback. By displaying all relevant process data on-screen, the laboratory workflow can be perfectly timed.








### Features

- Smart user interfaces
- Workflow oriented wizards
- Predefined templates
- Easy adjustment
- Quick start



### Place Load Teach Start

Multifaceted specimen carriers are available for various applications; Leica Microsystems' scanning templates are easily customized.



Mouse diaphragm muscle stained against neuro-filament 150. Mosaic: xyz: 5 x 5 x 101 images  
Green: Alexa Fluor 488-secondary antibody and acetylcholine receptors (Red: Alexa Fluor 647-alpha-bungarotoxin). Courtesy of Dr. R. Rudolf, Cellular Signaling in Skeletal Muscle, Karlsruhe Institute of Technology, Germany.

# Automated High Content Screening Simplifies D

LAS AF software simplifies routines. Leica Microsystems' goal is to make daily work as easy as possible so researchers can concentrate on the results, not on the imaging process.

To draw the correct conclusion, meta data is very important. LAS AF MATRIX M3 links meta data with the image to provide comprehensive information and enables researchers to go back to the single imaging results of each well at any time. The Leica HCS A data model ensures that accurate results are generated from acquisition to future image analysis.

## We Keep it Simple!

### **LAS AF MATRIX Mosaic Applications**

Fine details are as important as an information overview when evaluating experimental data. Leica HCS A includes re-engineered mosaic algorithms for excellent results at the push of a button.

Leica LAS AF MATRIX Mosaic automatically generates large high content images, providing both an overview and high resolution image at the same time.

Mouse diaphragm muscle stained against neuro-filament 150. Mosaic: xyz: 5 x 5 x 101 images. (Green: secondary antibody coupled to Alexa Fluor 488) and acetylcholine receptors (Red: alpha-bungarotoxin coupled to AlexaFluor 647). Courtesy of Dr. Rüdiger Rudolf, Cellular Signaling in Skeletal Muscle, Karlsruhe Institute of Technology, Germany.

## Get the Content!

### **LAS AF MATRIX Multiwell Applications**


Answers to more complex questions require more freedom in experimental design. Leica HCS A supports frequently used multiwell plate formats to automatically study multi-dimensional experiments. Time resolved or concentration dependent tests unveil true biological context.

Zebrafish, *Danio rerio*, Neurogenin - GFP. H2A  
Courtesy of J. Legradi, Dr. U. Liebel, KIT Karlsruhe Institute of Technology, Germany.




# aily Routines


High Resolution Single Image




High Content Mosaic



Fast Multiwell Plate Screening



High Content Information



A wide range of experiments are possible, with flexible screening conditions that can take care of even the smallest field scans.


# Gain Flexibility

## Adjustable Scanning Templates


You don't need to be a programmer – new scanning templates can be designed at micron scale to easily fit chamber slides, multi-well plates or simple spotter arrays. Once defined, the templates are ready to use for all applications and can even be conveniently shared between laboratories or communities.

## Imaging without Limits – MultiJob and MultiPositioning

Feel free to combine a variety of individual scan jobs for any area of interest of the specimen. The MultiJob – MultiPosition function provides maximum flexibility for your experiments. Several jobs, such as low resolution prescreens or multicolor 3D acquisition, can be freely combined. Change objective magnification or zoom in and out automatically. Individual settings can be adjusted for each position. From basic routines up to the most complex experiments, Leica HCS A greatly extends the spectra of applications.



Make it fit: adjustment flexibility for scanning templates.



Automated color coding assigned to the wells above provides easy control of multiple scan jobs.

## Trigger In – Automated Control

This function allows screening jobs immediately to start on external trigger signals. *Trigger In* provides new opportunities for external events to control internal process steps.

## Water Immersion Objective Control – Never Run Dry

A sufficient supply of immersion fluid is highly important for long-term observations. Water dispense volume, timing, and position can be controlled to maintain excellent optical conditions based on environmental conditions.



Automated water immersion objective

## Autofocus Routines

Several autofocus algorithms are available, optimized for different setups. The suitable routine is selected from a pull-down menu. After the initial screen, the software automatically creates a focus map with true sample topology. This map is used for fast, accurate z-positioning during the screening experiment. According to the size and planarity of the samples, the optimal number and positions of the autofocus points can be freely defined.

## Z-Drift Compensation

Live specimens can move or change size or shape over time, changing the z-position of interest. Microscope conditions can change due to temperature shifts. The algorithm adjusts the focus independently over time and provides sharp images throughout the experiment.



## Tracking Algorithm

As live organisms may change their xy-position, the center of intensity is calculated at each time point. If the single target is moving, the software automatically repositions the object of interest to the center of the objective, providing the best imaging conditions.

## Review On-the-Fly


Data is stored at predefined locations on a local hard disk or network storage device (NAS) via a TCP/IP protocol. Experiment data flows into a ring buffer to ensure that an unlimited stream of images enters the specified target folder.

The advantage: data analysis or review of image data is performed immediately. Image analysis starts when the experiment starts and provides fast results. Fast feedback loops automatically communicate with the microscope as it screens to modify the scanning parameters. This enables researchers to detect rare events when they happen.



### Autofocus Procedure

Z-stack images are acquired at freely selectable positions. The focus positions are determined and stored in a color coded focus map.



LAS AF live image viewer

# Automated High Content Screening for Quantification

## Detect Rare Events On-the-Fly!

### LAS AF MATRIX Developer Suite

Highly complex assay designs are now possible with LAS AF MATRIX Developer Suite. The Mitocheck project (1), conducted at the EMBL in Heidelberg, is an excellent example of comprehensive and flexible automation: The self-acting process automates mitosis identification.

### Pre-Screen – Object Identification

At first, a fast, low resolution pre-screen of sample plates is performed to identify the events of interest (1). After each scan, OME-image data are streamed on-the-fly to a buffer to be distributed online. The image information is stored on a network attached storage device (NAS) or server to be processed with the researcher's image analysis (2).

Fully automated, the pre-screen data is segmented, object features are extracted, and the target cells are classified. The object is selected and can clearly be identified by the meta data of the microscope (3). The position of interest is reported via a network protocol to the Computer Aided Microscopy (CAM) interface, which now automatically starts additional imaging protocols (4).

### High Content – Object Scan


In the second phase, high resolution scans start to acquire time resolved behavior of the target cells at the reported positions. The imaging mode automatically switches from pre-screen to high resolution acquisition. Now, multiple z-positions and additional channels are applied. Only the positions marked during the pre-screen are scanned. After each cycle, a pre-screen is repeated to find new cells entering mitosis, and the process starts again.

The method is highly efficient: all mitotic events can be identified and scanning of non-target cells is avoided. Toxic bleaching of cells that are not in mitosis is minimized. Thus, scanning speed is maximized at reduced laser exposure.

Microtubule secondary screen: scrambled sRNA Tubulin (green), H2B (red).  
Leica TCS SP5. Objective: 63x oil (pre-scan); zoomed, maximum projection:  
30 x 0.4  $\mu\text{m}$  slices, 2 channels (high resolution).  
Courtesy of Christian Conrad, EMBL, Heidelberg, Germany.  
(1) Mitocheck Project: [www.mitocheck.org](http://www.mitocheck.org).

## High Content Screening with

### Fast Pre-Screen




2



### Data Online

# Quantitative Assays on Leica TCS SP5

with Interactive System Control: Fully Automated Mitosis Acquisition



# The Perfect Match

## Data Interfaces




Open architecture for truly platform independent information exchange in an interactive environment – this is the goal attained by the new Leica HCS A data model.

### Experiment Meta Data Administration

Experiment IDs, description and meta data can be entered manually or by barcode. Additional experiment information can be added to the existing XML meta data file by external programming to provide comprehensive result data.

### Platform Independent Results Distribution


Leica HCS A imaging formats are platform independent and can be used on Apple MAC™ OS, Microsoft Windows® or LINUX1 platforms. The new Data Exporter provides OME-TIFF image files automatically, containing binary image data plus XML meta data structure.




Automated image analysis of multiple assay by DEFINIENS Cellenger®, Courtesy of Dr. R. Pepperkok, European Molecular Biology Laboratory (EMBL), Heidelberg, Germany.

## Image Analysis

Leica HCS A export formats are easily imported into all modern image analysis solutions such as e.g. DEFINIENS Cellenger®, ImageJ2, or MetaMorph®. Ensuring full compatibility to modern analysis platforms provides new options for target recognition, analysis and decision-making. Researchers benefit from applying existing algorithms or may create new ones to analyze data automatically.



Modern analysis tools offer solutions to identify, count, and measure target cells to obtain statistically relevant results. In addition, DEFINIENS Cellenger® quantifies relationships between target objects even in 2D or 3D. End users are free to choose among local, server-based, or clustered node analysis methods to maximize throughput and efficiently achieve high content screening results.



# Transparency of High Content Data


## Data Model

In the past, the original imaging and meta data had to migrate through a myriad of different conversion formats before ending up in a condensed Excel or Word document. Loss of data due to conversion is a problem of the past. The Leica export format follows the conventions of well-defined and well-formed structures, and can be read by all modern software platforms. Data conversion is no longer necessary. Transformation problems, data mix-up or transcription errors are avoided and processing time is saved. Additionally, even meta data can be combined with the results of Leica HCS A using external programs.


With Leica LAS AF MATRIX M3, researchers now get faster results because a clear picture of the experiment and result data is provided. From sample preparation to acquisition parameters to future image analysis data, you will never lose information with this scalable data model.


### Easy Collaboration – The Power of Networking

Sharing data is the key to successful collaborations. Leica HCS A realises its full power in a network environment. Information is exchanged in a LAN via basic TCP/IP protocols from the imaging system to network attached storage and analysis tools. Collaborative sharing of scientific data between laboratories in your facility provides an extendable, time saving way to achieve fast results.



### Data Model:





Fibroblasts Nuclei (DAPI, blue) and Actin (Phalloidin-TRITC red). Courtesy of Dr. G. Giese, MPI for Medical Research, Heidelberg, Germany.

**Annotations:**

MAC™ OS X is a registered trademark of Apple® Inc. Windows® is a registered trademark of the Microsoft® Corporation. (1) Linux is a free Unix-type operating system originally created by L. Torvalds with the assistance of developers around the world. Definiens® is a Registered Trademark of Definiens AG.

(2) ImageJ is a public domain Java image processing program inspired by National Institutes of Health, NIH Image for Windows®, Mac™ OS, Mac™ OS X and Linux. MetaMorph® is a Registered Trademark of MDS Analytical Technologies. Huygens Professional® is a Registered Trademark of SVI Scientific Volume Imaging.

(3) Open Microscopy Environment (OME) is a multi-site collaborative effort among academic laboratories and a number of commercial entities that produces open tools to support data management for biological light microscopy. Designed to interact with existing commercial software, all OME formats and software are free, and all OME source code is available under GNU public copyleft licenses. OME is developed as a joint project between research-active laboratories at the Dundee, NIA Baltimore, and Harvard Medical School and LOCI. LabVIEW™ is a registered trademark of NI National Instruments Inc. MATLAB™ is a registered trademark of The MathWorks™, Inc. Java™ is a registered trademarks of Sun Microsystems, Inc. C++ is a programming language standardized by ISO. C# is a programming language developed by Microsoft, Inc.


**The Benefit of Integration – Excellent Results**

Leica HCS A provides full connectivity within every laboratory environment by overcoming format barriers. Using existing software solutions saves time and money. Implementation is fast and ensures excellent research results for high content experiments.

# Automation Control via LAN


**Computer Aided Microscopy – Customize Your Imaging System**

Instant interaction, based on clear decisions is the key for success in science. Computer Aided Microscopy (CAM) provides the tool to immediately control a microscope system by full automation. A few easy-to-learn commands are sufficient to get you started.




**Get the Power**

The new CAM interface offers remote control of confocal and widefield systems, such as the Leica TCS SP5, TCS SPE, TCS LSI, Leica AF7000, AF6500 and AF6000 by LabVIEW, MATLAB or script based programming languages. Individual imaging jobs are quickly started, based on the decisions of image analysis, external trigger events or time loops.



Immediately after image acquisition, the data streams into the analysis tools for processing. A moment later, target cells are clearly classified and marked by the spatial position. Following the program, the instrument may now start a zoom-in or high resolution screen for more detailed observation. Due to the high speed of the process, even rare events are captured that previously may have been lost.





Overcome the limits of static experiments and immediately start individual imaging processes based on the results of image analysis. Unbiased algorithms perform objective target selections and automatically perform statistically relevant screens. Leica HCS A provides all the tools needed for excellent high content screening results.

### System Administration

The system follows standard rules of LAN administration to best fit into the facility's IT structure. The local administrator is always in full control of the system as LAS AF MATRIX M3 uses assigned permissions.

Standard remote system tests can be started by any user to ensure uptime and reproducible results. Additionally, Leica Microsystems offers technical service support by RemoteCare® instrument diagnosis via a secure internet protocol to minimize downtime.

# Powerful Platforms

## Hardware for High Content Screening

Leica Microsystems, world renowned for excellent high resolution image quality, offers confocal and widefield system platforms for automated high content screening. Choose the best tool for your applications.

### Leica TCS SP5

The only broadband confocal is a universal high-speed platform for parallel multi-channel micro imaging. Leica Microsystems' AOBs (Acousto-Optical Beam Splitter) technology provides full spectral detection and highest transmission.

### Leica TCS SPE


The compact, robust confocal system is cost-effective and extremely easy to use. The system performs sequential multi-channel micro image acquisition analysis up to eight colors. A special glass prism and precise spectral selectors provide maximum spectral efficiency.

### Leica TCS LSI


The platform provides high resolution imaging from micro to macro. The combination of confocal and optical zoom offers maximum flexibility for specimens up to 16 mm. Sequential multi-channel micro plus macro image acquisition analysis up to eight colors. A special glass prism provides the transmission for spectral analysis.

### Leica AF7000, AF6500 and AF6000


The Leica widefield systems AF7000, AF6500 and AF6000 offer a large portfolio of different CCD and EM CCD cameras for fastest acquisition speed and highest quantum efficiency. Different widefield illumination options are possible. The widefield screening systems are based on the inverted microscope DMI6000 B or on the upright microscope DM6000 B. An ideal system solution for primary and secondary screens.



Drosophila, Leica TCS SP5



Fibroblasts, Leica TCS SPE



Mouse embryo, Leica TCS LSI



Leica TCS SP5



Leica TCS SPE


# Gain Scientific Advantage

## Intelligent Automation

The Leica HCS A technology provides an excellent platform for fast, efficient high content screening. Leica combines highest resolution for all specimen sizes with LAS AF MATRIX M3 automation software. Maximum application flexibility with varying levels of automation creates unrivaled freedom of experimental design for today and tomorrow.

Leica HCS A can interface to every laboratory environment in an optimal way via standardized interfaces. Open, well-defined architecture and Open Microscopy Environment (OME) formats are compatible with many image analysis technologies. The perfect match between scalable export formats and platform independent device control via Computer Aided Microscopy (CAM) interfaces creates added value. Obtain more results in less time – efficiently and statistically verified.

Gain scientific advantage from Leica HCS A for high content screening and amplify the power of imaging.



HeLa Cells expressing Histone H2B and Tubulin. Cells observed on Leica AF7000 using the Water Immersion Micro Dispenser and HCX PL APO 63x/1.20 W CORR CS objective over 3 days in 37°C. Courtesy of Dr. Daniel Gerlich (ETH Zurich) Jutta Bulkescher and Dr. Stefan Terjung (ALMF, EMBL Heidelberg).

### Acknowledgements

We gratefully acknowledge the scientists cited in the brochure and below for providing images:

#### 1, 4 *Drosophila melanogaster* (egg chamber)

Green: Actin, Alexa 488-Phalloidin; Red: Cortex, Egalitarian; Red Blue: hnRNP, Cy5; Grey: Nuclei, DAPI.

Courtesy of Sonja Lopez de Quinto, Florence Besse and Oliver Hachet, EMBL, Heidelberg, Germany.

#### 2, 6 Mouse embryo


Nuclear and neurofilament staining of mouse embryo (10.5 days post coitum). Courtesy of ICI, Imaging Center of IGBMC, Strasbourg, France.

#### 3 COS 7 cells

Green: uncharacterized protein, GFP; Red:  $\alpha$ -Tubulin, Cy3; Blue: Nuclei, DAPI. Courtesy of Prof. Wei Bian, Cell Research Center, Institute of Biochemistry and Cell Biology, SIBS, CAS, Shanghai, China.

#### 5 Mouse fibroblasts

Green: F-Actin, FITC; Red: Tubulin, Cy5; Blue: Nuclei, DAPI. Courtesy of Dr. Günter Giese, Max Planck Institute for Medical Research, Heidelberg, Germany.



Leica TCS LSI



Leica AF7000 based on DMI6000 B

# “With the user, for the user”

## Leica Microsystems

Leica Microsystems operates globally in four divisions, where we rank with the market leaders.

### • Life Science Division

The Leica Microsystems Life Science Division supports the imaging needs of the scientific community with advanced innovation and technical expertise for the visualization, measurement, and analysis of microstructures. Our strong focus on understanding scientific applications puts Leica Microsystems' customers at the leading edge of science.

### • Industry Division

The Leica Microsystems Industry Division's focus is to support customers' pursuit of the highest quality end result. Leica Microsystems provide the best and most innovative imaging systems to see, measure, and analyze the microstructures in routine and research industrial applications, materials science, quality control, forensic science investigation, and educational applications.

### • Biosystems Division

The Leica Microsystems Biosystems Division brings histopathology labs and researchers the highest-quality, most comprehensive product range. From patient to pathologist, the range includes the ideal product for each histology step and high-productivity workflow solutions for the entire lab. With complete histology systems featuring innovative automation and Novocastra™ reagents, Leica Microsystems creates better patient care through rapid turnaround, diagnostic confidence, and close customer collaboration.

### • Surgical Division

The Leica Microsystems Surgical Division's focus is to partner with and support surgeons and their care of patients with the highest-quality, most innovative surgical microscope technology today and into the future.

The statement by Ernst Leitz in 1907, “with the user, for the user,” describes the fruitful collaboration with end users and driving force of innovation at Leica Microsystems. We have developed five brand values to live up to this tradition: Pioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement. For us, living up to these values means: **Living up to Life.**

### Active worldwide

Australia:	North Ryde	Tel. +61 2 8870 3500	Fax +61 2 9878 1055
Austria:	Vienna	Tel. +43 1 486 80 50 0	Fax +43 1 486 80 50 30
Belgium:	Groot Bijgaarden	Tel. +32 2 790 98 50	Fax +32 2 790 98 68
Canada:	Richmond Hill/Ontario	Tel. +1 905 762 2000	Fax +1 905 762 8937
Denmark:	Herlev	Tel. +45 4454 0101	Fax +45 4454 0111
France:	Rueil-Malmaison	Tel. +33 1 47 32 85 85	Fax +33 1 47 32 85 86
Germany:	Wetzlar	Tel. +49 64 41 29 40 00	Fax +49 64 41 29 41 55
Italy:	Milan	Tel. +39 02 574 861	Fax +39 02 574 03392
Japan:	Tokyo	Tel. +81 3 5421 2800	Fax +81 3 5421 2896
Korea:	Seoul	Tel. +82 2 514 65 43	Fax +82 2 514 65 48
Netherlands:	Rijswijk	Tel. +31 70 4132 100	Fax +31 70 4132 109
People's Rep. of China:	Hong Kong	Tel. +852 2564 6699	Fax +852 2564 4163
Portugal:	Lisbon	Tel. +351 21 388 9112	Fax +351 21 385 4668
Singapore		Tel. +65 6779 7823	Fax +65 6773 0628
Spain:	Barcelona	Tel. +34 93 494 95 30	Fax +34 93 494 95 32
Sweden:	Kista	Tel. +46 8 625 45 45	Fax +46 8 625 45 10
Switzerland:	Heerbrugg	Tel. +41 71 726 34 34	Fax +41 71 726 34 44
United Kingdom:	Milton Keynes	Tel. +44 1908 246 246	Fax +44 1908 609 992
USA:	Bannockburn/Illinois	Tel. +1 847 405 0123	Fax +1 847 405 0164

and representatives in more than 100 countries