

POWER UP ROUTINES WITH EFFICIENCY AND COMFORT

Visoria B laboratory microscope

The Visoria B laboratory microscope is for users in life science and clinical laboratories who perform applications to support human, animal and plant histology, and pathology examinations.

UPGRADE YOUR CURRENT MICROSCOPE

If you already have a DM2500, DM2700 M, or DM2700 P upright microscope, you can upgrade to a Visoria stand without paying the full price of an entirely new microscope. This means you will not need a new stage, objectives, filters, tubes or eyepieces - your current components are fully compatible with Visoria.

The Visoria stand is available in either a traditional (eyepieces) or digital configuration.

Your advantages with a Visoria stand:

- Save time with optimized light settings. The light management function automatically adjusts brightness when changing magnification or contrast methods.
- Operate your microscope with ease With color-coding quickly match the proper aperture diaphragm and objective, a built-in focus stop helps protect specimens and objectives, and a three-gear system enables finer focusing at higher magnification.
- Reduce the risk of discomfort and repetitive strain injury when spending long hours at the microscope, you can minimize repetitive movements with the Visoria series by adjusting the height and torque of the stage control.
- The knobs for controlling the XY stage position and focus of the Visoria microscope(s) are in close proximity and easily adjustable with one hand.
- The microscope can be easily changed between rightand left-handed operation. Especially beneficial when the microscope is shared with other users.
- Simplify your documentation. Image capture with the press of a button and system settings automatically stored with images.



Powered by Enersight software

The Enersight software platform helps you compare, measure, and share seamlessly with a single intuitive interface.

- Observe samples with a larger field of view and higher resolution using the 2D Stitching function.
- Acquire sharp images of samples with extended depth of field (EDOF).
- Choose your optimal illumination and camera parameters with the Best Image function.



Specifications Visoria B

Microscope	
Size and weight	Length: 410 mm, 456 mm (with fluorescence lamphouse), width: 331 mm, height: 505 mm, approx. 18 kg
	(depending on configuration)
Stand	Illumination toggle buttons with status indicator, image capture button, antimicrobial surface with
	AgTreat according to ISO 22196
Optics	
Nosepiece	Encoded 6x (M25)
Eyepieces (FOV)	20 / 22 / 25 mm
Tubes	Wide range of standard, ergonomic and phototubes with different beam splitters available
Digital version	Digital version with 10" screen / tablet
Ergonomic accessories	Wide range of ergonomic accessories available (ErgoTubes, ErgoLift, ErgoModules)
Encoded illumination	IL and TL: High-power white LED, encoded 4-color fluorescence illumination: emission wavelengths:
management	UV 370 / 12 nm, B 465 / 20 nm, G 546 / 82 nm, R 624 / 15 nm, further external light sources on request
	(non-encoded)
Fluorescence light axis	Optional
Incident light (IL)	Methods: Brightfield, fluorescence, DIC, qualitative polarization
Transmitted light axis	Manual, fixed and flip-top condenser operation with color-coded diaphragm assistant
Transmitted light (TL)	Methods: Brightfield, darkfield, phase contrast, DIC, qualitative polarization
Operation	
Stage	Manual XY stage with ceramic plate, additional manual xy stages available, 1-specimen or 2-specimen holder
Stage control	Left-, right-handed stage, torque-adjustable handle
Focus drive	Height-adjustable focus knobs, 19 mm travel range, maximum 28 mm total stage stroke depending on
	stage and condenser type, 2-gear focus drive (coarse / fine) with 1 mm scale, 3-gear focus drive
	(coarse / medium / fine) with 140, 4 and 1 µm scale, torque adjustment, and adjustable upper focus-sto
Accessories	
Analyzer	Fixed, 180°, 360°
Polarizer	Fixed and rotatable
General specifications	
Supply voltage	100–240 V AC, 50 / 60 Hz, power consumption max. 15 W
Ambient conditions	15–35°C, relative humidity max. 80% up to 30°C (non-condensing)

CONNECT WITH US!



Leica Microsystems GmbH | Ernst-Leitz-Strasse 17–37 | D-35578 Wetzlar (Germany) Tel. +49 (0) 6441 29-0 | F +49 (0) 6441 29-2599

www.leica-microsystems.com

